# ELECTRICITY/ELECTRONICS

The introduction to *Electricity/Electronics* for the students usually begins with experiments such as the simple electric circuit. They learn the basic principles that are essential for understanding the experiments based on them.

The use of the innovative, modular plug-in board makes it particularly easy to reduce the experiments to a minimum. The safety sockets on the plug-in board enable the use of safety connecting leads. *Five* Electrics Sets provide six topic areas with 154 experiments.

The measurements of current and voltage can be carried out both with the Mobile-CASSY 2 WiFi and with multimeters.



## LP3.1.7.1C Equipotential lines between identically shaped electrodes

In this experiment, students record the equipotential lines between two identically shaped electrodes by searching for the points of equal potential difference between the 0 V electrode and the steel needle with the Mobile-CASSY 2 WiFi and depict these points on a graph. For this experiment you will need the set Science Lab Electrics EL1 (207 131S).

Further information about our curriculum-compliant topics and student experiments as well as the corresponding sets can be found on the following pages.

# **Overview** of topics and sets

EXPERIN	IENT TOPICS	REQUIRED SETS	NO. EXPERIMENTS	DETAILS
LP3.1	ELECTROSTATICS			
LP3.1.1	CONTACT ELECTRICITY	Electrics EL1	25	PAGE 72
LP3.1.2	FORCES ACTING BETWEEN CHARGES			
LP3.1.3	ELECTROSTATIC INDUCTION			
LP3.1.4	CHARGE STORAGE			
LP3.1.5	ELECTROSTATIC INTERACTION			
LP3.1.6	INSULATORS AND CONDUCTORS	207 131S		
LP3.1.7	EQUIPOTENTIAL LINES			
LP3.1.8	PLATE CAPACITOR			
LP3.2	MAGNETISM			
LP3.2.1	MAGNETIC FORCES AND FIELDS	Electrics EL2	12	PAGE 78
LP3.2.2	MAGNETIC INDUCTION			
LP3.2.3	MAGNETIC FIELDS			
		207 132S		
LP3.3	BASIC ELECTRICAL CIRCUITS AND ELECTROCH	EMISTRY		
LP3.3.1	ELECTRICAL CIRCUITS AND SWITCHES	Electrics EL3	40	PAGE 84
LP3.3.2	ELECTRICAL MEASUREMENT METHODS			
LP3.3.3	OHMIC RESISTANCE			
LP3.3.4	SPECIAL RESISTORS			
LP3.3.5	VOLTAGE SOURCES			
LP3.3.6	ELECTRICAL APPLICATION CIRCUITS	207 133S		
LP3.3.7	ELECTROCHEMISTRY			
LP3.4	ELECTROMAGNETISM AND INDUCTION	,		
LP3.4.1	ELECTROMAGNETISM	Electrics EL3 Electrics EL4	21	PAGE 90
LP3.4.2	ELECTROMAGNETIC APPLICATIONS			
LP3.4.3	INDUCTION			-
LP3.4.4	TRANSFORMERS			-
LP3.4.5	APPLICATIONS OF INDUCTION			-
LP3.4.6	COILS IN DIRECT AND ALTERNATING CURRENT CIRUITS	207 1335 207 134S		
LP3.5	MOTORS AND GENERATORS			
LP3.5.1	GENERATORS	Electrics EL3 Electrics EL4 Electrics EL5	14	PAGE 96
LP3.5.2	ELECTRIC MOTORS			
		207 133S 207 134S 207 135S		
EXPERIM		REQUIRED SETS	NO EXPERIMENTS	DFTAILS

S	
$\overline{=}$	
D	
R	
5	
щ	

		207 1335	207 1345	207 1355		
EXPERIM	IENT TOPICS		REQUIRED SET	S	NO. EXPERIMENTS	DETAILS
LP4.1	BASIC ELECTRONIC CIRCUITS					
LP4.1.1	CAPACITORS	Electrics E	EL3 I	Electrics EL6	42	PAGE 102
LP4.1.2	RELAY CIRCUITS					
LP4.1.3	DIODES		Т			
LP4.1.4	TRANSISTORS		/ T (	-		
LP4.1.5	DIODE CIRCUITS					
LP4.1.6	FLIP-FLOPS	207 133	S	207 136S		
LP4.1.7	AMPLIFIER CIRCUITS					
LP4.1.8	SOLAR CELLS					

## ELECTRICS – EL1

PHYSICS

### OVERVIEW OF OUR CURRICULUM-COMPLIANT EXPERIMENTS

### LP3.1 ELECTROSTATICS

LP3.1.1 Contact electricity	
LP3.1.1.1Proof of charge types on friction rods with a glow lampLP3.1.1.2Proof of charge types on foils and sheets with a glow lampLP3.1.1.3Contact electricity generated by friction	
LP3.1.2 Forces acting between charges	
LP3.1.2.1Forces acting between charges on friction rodsLP3.1.2.2Forces acting on a charged pendulum pairLP3.1.2.3Operating principle of an electroscope	
LP3.1.3 Electrostatic induction	
LP3.1.3.1Induction phenomena with conductors and non-conductorsLP3.1.3.2Induction phenomena of a pointerLP3.1.3.3Electric induction on a pair of pendulumsLP3.1.3.4Induction phenomena on an electroscope caused by friction rods	
LP3.1.4 Charge storage	
LP3.1.4.1Conductive bodies as charge storage devicesLP3.1.4.2Proof of charges on a Faraday cup	
LP3.1.5 Electrostatic interaction	
LP3.1.5.1Electrostatic forces between friction rod and pendulumLP3.1.5.2Charge transfer through a pendulum	
LP3.1.6 Insulators and conductors	
LP3.1.6.1Charges on insulatorsLP3.1.6.2Proof of conductivity with an electroscopeLP3.1.6.3Point discharge	
LP3.1.7 Equipotential lines	
LP3.1.7.1Equipotential lines between identically shaped electrodesLP3.1.7.1CEquipotential lines between identically shaped electrodes (with Mobile-CASSY 2 WiFi)LP3.1.7.2Equipotential lines between non-identically shaped electrodesLP3.1.7.2CEquipotential lines between non-identically shaped electrodes (with Mobile-CASSY 2 WiFi)LP3.1.7.3Distortion of equipotential linesLP3.1.7.4Distortion of equipotential lines	OIGITAL OIGITAL
19.1.7.50     Distriction of equipotential lines (with Mobile-CASSE2 with)       1P3.1.8     Plate capacitor	DIGITAL
LP3.1.8.1 The electric field in a plate capacitor	
LP3.1.8.1C The electric field in a plate capacitor (with Mobile-CASSY 2 WiFi)	DIGITAL

EXPERIMENTS



## LP3.1.4.2 Proof of charges on a Faraday cup

Students will show that the two different charges can be retained on a Faraday cup and detected with the glow lamp. For this experiment you will need the set **Science Lab Electrics EL1 (207 131S)**.

### OVERVIEW OF EQUIPMENT REQUIRED FOR PERFORMING EXPERIMENTS





### Science Lab Electrics EL1 (Set)

Student experiment set of the student experiment system Science Lab in the field of physics. Set-up material for one working group in pre-formed tray. With the equipment set EL1 (207 131S), 25 experiments at school, college and university level for worldwide curriculums can be performed. The students deal with the topics electrostatics and electric fields. While working out the curriculum required topics, they are also trained in communication and assessment skills. In combination with the Mobile-CASSY 2 WiFi (524 005W), there are additional evaluation options which enable the students digital learning.

### Scope of delivery:

Count	Name	Count	Name
2	Safety connecting lead 50 cm, red	1	Clamping plug
2	Safety connecting lead 50 cm, blue	1	Tray, high
3	Adapter 4-mm plug/4-mm socket	1	Lid for tray
1	Crocodile clip, polished	1	Cord
1	Glow lamp, tubular 90 V	1	Plastic clips, pair
1	Electroscope	1	Electrostatic pendulums, pair
2	Friction rods, PVC and acrylic	1	Microfibre cloth
1	Induction plate 8 cm x 4 cm	1	Steel needle
1	Bar electrodes for 54509, set 2	5 out of	Acetate foils 300 x 300 x 0.1 mm, set of 10
1	Set of 3 round electrodes for 54509	207 1315	Science Lab Electrics EL1 (Set)
1	Faraday's cup		

### ADDITIONALLY REQUIRED TO PERFORM ALL EXPERIMENTS

Count	CatNo.	Name	Description
1	521 487	AC/DC Power supply PRO 012 V/3 A	
1	524 005W	Mobile-CASSY 2 WiFi	for digital experiments
1	531 120	Multimeter LDanalog 20	alternative for analog measurements
1 \dditi	531 120 onally req	Multimeter LDanalog 20 uired per <mark>class</mark>	alternative for analog measurements
1 <b>\dditi</b> Count	531 120 onally req CatNo.	Multimeter LDanalog 20 uired per <mark>class</mark> Name	alternative for analog measurements Description
1 Additi Count 1	531 120 onally req CatNo. 520 713	Multimeter LDanalog 20 uired per class Name LIT: LP3 Science Lab Electricity, digital	alternative for analog measurements Description



# **OVERVIEW OF ADVANTAGES**

- Students build their "own" electroscope and learn about its function
- Straightforward experiments for the detection and storage of different charges and electric fields
- All electric and electronic components are in transparent housings for a safe and long term use and labelled with the same electronic symbols as real life devices
- Acquired skills: Design and function of electrical measuring instruments

### STUDENT MEASURING DEVICE

DIGITAL CLASS / EDUCATION



Mobile-CASSY 2 WiFi

The universal student measuring device with WiFi for all measuring tasks in physics, chemistry and biology.

524 005W Mobile-CASSY 2 WiFi

You can find detailed information on the Mobile-CASSY 2 WiFi on page 228.







With the Mobile-CASSY 2 WiFi, voltage (U), current (I), power (P) and energy (E) can be measured via 4 mm safety sockets.

### LITERATURE PACKAGES

Here you will find an overview of our literature packages.

You can find detailed information on our literature on the internet at www.leybold-shop.com.









### Printed version available in ring file LIT: LP3.1 Electrostatics

Detailed experiment instructions relating to Science Lab Set EL1 (207 131S). Describes 25 experiments from the field of electrostatics.

Topics:

Contact electricity; Force acting between charges; Electrostatic induction; Charge storage; Electrostatic interaction; Insulators and conductors; Equipotential lines; Plate capacitor

520 7131EN

LIT: LP3.1 Electrostatics

# LIT: LP3 Science Lab Electricity, digital includes only ONE subject area

Comprehensive physics experiment instructions in the field of electricity for the Science Lab. Contains 154 experiments on electrostatics, magnetism, basic electrical circuits and electrochemistry, electromagnetism and induction, motors and generators and electronics.

Includes all interactive experiment instructions (Lab Docs) as html file.

LIT: LP3 Science Lab Electricity, digital 520 713

### LIT: LP Science Lab Physics, digital

Comprehensive physics experiment instructions for the Science Lab. Contains 450 experiments in the fields of mechanics, energy, electricity and electronics, optics, atomic and nuclear physics.

includes ALL subject areas

Includes all interactive experiment instructions (Lab Docs) as html file.

#### 520 71 LIT: LP Science Lab Physics, digital

### Technical data of the digital version:

- Product key for literature (activation & selection of one literature language in LeyLab)
- Can then be used in LeyLab and Document Center (school/institute licence)
- System requirements: Document Center:

- PC with Windows 7 or higher; internet access during installation; local network for distribution to students LevLab:

- PC, tablet or smartphone with a current browser; internet access

### ADDITIONAL STORAGE ACCESSORIES



You can find detailed information on additional storage accessories from page 228.

www.ld-didactic.com

### INTRODUCING THE TOPIC

## Experiment and research – Discover electrostatic phenomena with fun

Even though the names of the experiments may not sound like fun to students first, the electrostatics experiments vividly illustrate many everyday life experiences.







## SOME EXAMPLES OF STUDENT QUESTIONS THAT WILL BE ANSWERED:

- Why do I get an electric shock if I touch a metal door after I have walked on carpet with shoes on?
- Why does that seem to be worse with certain shoes than with others?
- Why does hair stick to a balloon if I rub the balloon on my t-shirt first?

# ELECTRICS – EL2

PHYSICS

### OVERVIEW OF OUR CURRICULUM-COMPLIANT EXPERIMENTS

## LP3.2 MAGNETISM

LP3.2.1	Magnetic forces and fields
LP3.2.1.1 LP3.2.1.2 LP3.2.1.3	Magnetic and non-magnetic materials Position of the magnetic poles on a bar magnet Polarity of magnets
LP3.2.2	Magnetic induction
LP3.2.2.1 LP3.2.2.2 LP3.2.2.3	Magnetisation Disassembling magnets Combining magnets
LP3.2.3	Magnetic fields
LP3.2.3.1 LP3.2.3.2 LP3.2.3.3 LP3.2.3.4 LP3.2.3.5 LP3.2.3.6	Demonstration of a magnetic field with iron powder Field lines of a bar magnet Model experiment on the earth's magnetic field Field lines of a horseshoe magnet Field lines of attracting magnetic poles Field lines of repelling magnetic poles



LP3.2.1.1 Magnetic and non-magnetic materials



LP3.2.3.1 Demonstration of a magnetic field with iron powder

### OVERVIEW OF EQUIPMENT REQUIRED FOR PERFORMING EXPERIMENTS







### Science Lab Electrics EL2 (Set)

Student experiment set of the student experiment system Science Lab in the field of physics. Set-up material for one working group in pre-formed tray. With the equipment set EL2, 12 experiments at school, college and university level for worldwide curriculums can be performed. The students deal with the topic magnetism. While working out the curriculum required topics, they are also trained in communication and assessment skills. In combination with the Mobile-CASSY 2 WiFi (524 005W), there are additional evaluation options which enable the students digital learning.

### Scope of delivery:

Count	Name	Count	Name
2	Bar magnet	1	Magnetic field indicator
2	Plotting compass	1	Shaker for iron filings
1	Magnetizable rods, set of 4	1	Tray, Iow
1	Hemisphere for earth's magnetism	1	Iron powder, reduced, 50 g
1	Magnetism experimental field	1	Cord
1	Iron yokes, pair	207 1325	Science Lab Electrics EL2 (Set)
1	Pocket compass		

### ADDITIONALLY REQUIRED TO PERFORM ALL EXPERIMENTS

Count	onally requ CatNo.	ired per working group Name	Description
1	510 55	Direction-finding compass	
Additio	onally requ	lired per class	
Count	CatNo.	Name	Description
1	520 713	LIT: LP3 Science Lab Electricity, digital	



DIGITAL CLASS / EDUCATION

# **OVERVIEW OF ADVANTAGES**

- Understandable experiments to distinguish between attracting and repelling magnetic forces which can be visualised with the help of the magnetism experimental field and iron powder
- Investigation of the Earth's magnetic field in a simple model
- Acquired skills: Explanation of everyday phenomena using simple models

### STUDENT MEASURING DEVICE



### Mobile-CASSY 2 WiFi

The universal student measuring device with WiFi for all measuring tasks in physics, chemistry and biology.

524 005W Mobile-CASSY 2 WiFi

You can find detailed information on the Mobile-CASSY 2 WiFi on page 228.



**LEYBOLD**<sup>®</sup>





You can find detailed information on additional storage accessories from page 228.

### INTRODUCING THE TOPIC

# Discover and experience magnetism by yourself



## DESCRIBING FIELD LINES AND INTERPRETING THE BEHAVIOUR OF MAGNETS





## MAKING MAGNETIC FIELDS VISIBLE



DRAWING CONCLUSIONS ABOUT THE EARTH'S MAGNETIC FIELD



## UNDERSTANDING THE FUNCTIONS OF A COMPASS

USING THE HEMISPHERE TO UNDERSTAND THAT THE EARTH IS A DIPOLE

# ELECTRICS – EL3

OVERVIEW OF OUR CURRICULUM-COMPLIANT EXPERIMENTS

### LP3.3 BASIC ELECTRICAL CIRCUITS AND ELECTROCHEMISTRY

LP3.3.1	Electrical circuits and switches		
LP3.3.1.1 LP3.3.1.2 LP3.3.1.3 LP3.3.1.4 LP3.3.1.5	The simple circuit Conductors and non-conductors (insulators) Switching over Two-way switches AND gate, OR gate		
LP3.3.2	Electrical measurement methods		
LP3.3.2.1 LP3.3.2.1C LP3.3.2.2 LP3.3.2.2C	Measuring current intensity in a simple circuit Measuring current intensity in a simple circuit (with Mobile-CASSY 2 WiFi) Measuring voltage in a simple circuit Measuring voltage in a simple circuit (with Mobile-CASSY 2 WiFi)		DIGITAL
LP3.3.3	Ohmic resistance		
LP3.3.3.1 LP3.3.3.1C LP3.3.3.2	Ohm's law Ohm's law (with Mobile-CASSY 2 WiFi) How a wire's resistance depends on its material, length and cross-section		DIGITAL
LP3.3.3.2C LP3.3.3.3 LP3.3.3.3C	How a wire's resistance depends on its material, length and cross-section (with Mobile-CASSY 2 WiFi) Voltage distribution in a current-carrying wire (potentiometer) Voltage distribution in a current-carrying wire (potentiometer) (with Mobile-CASSY 2 WiFi)		DIGITAL
LP3.3.3.4 LP3.3.3.4C LP3.3.3.5	Resistors in series Resistors in series (with Mobile-CASSY 2 WiFi) Resistors in parallel		DIGITAL
LP3.3.3.5C	Kesistors in parallel (with Mobile-CASSY 2 WiFi)		DIGITAL
	Special capacitors		
LP3.3.4.1 LP3.3.4.1C LP3.3.4.2 LP3.3.4.2C	Temperature-dependent resistors (NTC) Temperature-dependent resistors (NTC) (with Mobile-CASSY 2 WiFi) Light-dependent resistors LDR (photo-conductive cell) Light-dependent resistors LDR (photo-conductive cell) (with Mobile-CASSY 2 WiFi)		DIGITAL
LP3.3.5	Voltage sources		
LP3.3.5.1 LP3.3.5.1C LP3.3.5.2	Parallel and series connection of monocells Parallel and series connection of monocells (with Mobile-CASSY 2 WiFi) Terminal voltage and internal resistance of a voltage source		DIGITAL
LP3.3.5.2C	Terminal voltage and internal resistance of a voltage source (with Mobile-CASSY 2 WiFi)		DIGITAL
LP3.3.0	Electrical application circuits		
LP3.3.6.1C LP3.3.6.2 LP3.3.6.3 LP3.3.6.4	Self-heating and temperature sensitivity in wire-wound resistors (with Mobile-CASSY 2 WiFi) Model of a fuse Bimetal switch (model of a fire alarm) Power and work of an electrical current		DIGITAL
LP3.3.6.4C	Power and work of an electrical current (with Mobile-CASSY 2 WiFi)		DIGITAL
LP3.3.7	Electrochemistry		
LP3.3.7.1 LP3.3.7.1C	Conductivity of aqueous solutions (electrolytes) Conductivity of aqueous solutions (electrolytes) (with Mobile-CASSY 2 WiFi) Belation between current and voltage in an electrolyte		DIGITAL
LP3.3.7.2C	Relation between current and voltage in an electrolyte (with Mobile-CASSY 2 WiFi)		DIGITAL
LP3.3.7.3 LP3.3.7.4 LP3.3.7.4C	Galvanising Galvanic cells Galvanic cells (with Mobile-CASSY 2 WiFi)	40	DIGITAL
For experiments ma	rked with "C", the measurements are carried out digitally with the Mobile-CASSY 2 WiFi.	EXPERIMENTS	



LP3.3.1.4 Two-way switches

### OVERVIEW OF EQUIPMENT REQUIRED FOR PERFORMING EXPERIMENTS



### SCIENCE LAB PHYSICS - ELECTRICITY/ELECTRONICS

![](_page_16_Picture_1.jpeg)

![](_page_16_Picture_2.jpeg)

### Science Lab Electrics EL3 (Set)

Student experiment set of the student experiment system Science Lab in the field of physics. Set-up material for one working group in pre-formed tray. With the equipment set EL3, 40 experiments at school, college and university level for worldwide curriculums can be performed. The students deal with the topics electrical basic circuits and electrochemistry. While working out the curriculum required topics, they are also trained in communication and assessment skills. In combination with the Mobile-CASSY 2 WiFi (524 005W), there are additional evaluation options which enable the students digital learning.

Scone	of	deliverv	
JCOPC	01	uclively	

Count	Name
2	Coupling plug 4 mm
1	Bimetallic strip
2	Safety connecting lead 50 cm, red
2	Safety connecting lead 50 cm, blue
1	Safety connecting lead 50 cm, black
1	Bridging plugs STE 2/19, set of 10
1	Adapter 4-mm plug/4-mm socket
2	Crocodile clip, polished
1	Conductors/insulators, set of 6
1	Wrapping plate for wires
2	Plug-in board safety socket, 20/10
2	Monocell holder STE 2/50
1	Resistor 47 Ohm, STE 2/19
2	Resistor 100 Ohm, STE 2/19
1	Resistor 1 kOhm, STE 2/19
1	Resistor 10 kOhm, STE 2/19

Count	Name	
1	Variable resistor 47 kOhm, STE 2/19	
1	Photoresistor LDR 05, STE 2/19	
1	NTC resistor 2.2 kOhm, STE 2/19	
1	PTC resistor 100 Ohm, STE 2/19	
2	Lamp holder E10, lateral, STE 2/19	
1	Toggle switch STE 2/19	
2	Plug-in holder STE	
1	Contact strip	
2	Change-over switch STE 4/50	
1	Flat trough/electrolysis cell	
2	Plate electrode copper 76 x 40 mm	
1	Plate electrode zinc 76 x 40 mm	
1	Plate electrode iron 76 x 40 mm	
1	Tray, high	
1	Grindstone	
207 133S	Science Lab Electrics EL3 (Set)	

ADDITIONALLY REQUIRED TO PERFORM ALL EXPERIMENTS

![](_page_16_Picture_9.jpeg)

![](_page_16_Picture_10.jpeg)

# **OVERVIEW OF ADVANTAGES**

- The plug-in system enables even larger circuits to be mounted on the plug-in boards
- Variable plug-in board with safety sockets for 4 mm plugs
- Set up experiments in L- or T-shape with the plug-in board
- Easily expandable for more complex circuits by plugging several boards together
- Wrapping plate for wires allows easy experimentation on the resistance of wires while using less resources

### STUDENT MEASURING DEVICE

![](_page_17_Picture_8.jpeg)

![](_page_17_Picture_9.jpeg)

### Mobile-CASSY 2 WiFi

The universal student measuring device with WiFi for all measuring tasks in physics, chemistry and biology.

524 005W Mobile-CASSY 2 WiFi

You can find detailed information on the Mobile-CASSY 2 WiFi on page 228.

### ADDITIONALLY REQUIRED TO PERFORM ALL EXPERIMENTS

Auuiti	ionally req	uired per student	
Count	CatNo.	Name	Description
1	610 010	Laboratory safety goggles, Focomax	
Additi	ionally req	uired per working group	
Count	CatNo.	Name	Description
1	505 07	Bulbs, 4 V/0,16 W, E10, Set of 10	
1	505 08	Bulbs, 12 V/3 W, E10, Set of 10	
2	505 11	Bulbs, 2.5 V/0,25 W, E10, Set of 10	
1	521 487	AC/DC Power supply PRO 012 V/3 A	
1	524 005W	Mobile-CASSY 2 WiFi	for digital experiments
2	531 120	Multimeter LDanalog 20	alternative for analog measurements
2	685 48	Mono cell 1.5 V (IEC R20)	Voltage sources experiments (LP3.3.5)
A 1 11-1			
Additi	ionally req	uired per class	
Additi Count	ionally req CatNo.	uired per <mark>class</mark> Name	Description
Additi Count 1	ionally req CatNo. 520 713	uired per class Name LIT: LP3 Science Lab Electrics, digital	Description
Additi Count 1 1	ionally req CatNo. 520 713 550 42	uired per class Name LIT: LP3 Science Lab Electrics, digital Constantan resistance wire, 0.35 mm diameter, 100 m	Description Ohmic resistance experiments (LP3.3.3)
Additi Count 1 1 1	ionally req CatNo. 520 713 550 42 550 46	uired per class Name LIT: LP3 Science Lab Electrics, digital Constantan resistance wire, 0.35 mm diameter, 100 m Chrome-nickel resistance wire, 0.25 mm diameter, 100 m	Description         Ohmic resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)
Additi Count 1 1 1 1	CatNo. 520 713 550 42 550 46 550 47	uired per class Name LIT: LP3 Science Lab Electrics, digital Constantan resistance wire, 0.35 mm diameter, 100 m Chrome-nickel resistance wire, 0.25 mm diameter, 100 m Chrome-nickel resistance wire, 0.35 mm diameter, 100 m	Description         Ohmic resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)
Additi Count 1 1 1 1 1 1	CatNo. 520 713 550 42 550 46 550 47 550 51	uired per class Name LIT: LP3 Science Lab Electrics, digital Constantan resistance wire, 0.35 mm diameter, 100 m Chrome-nickel resistance wire, 0.35 mm diameter, 100 m Iron resistance wire, 0.35 mm diameter, 100 m	Description         Ohmic resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)
Additi Count 1 1 1 1 1 1 1	CatNo. 520 713 550 42 550 46 550 47 550 51 672 9650	uired per class Name LIT: LP3 Science Lab Electrics, digital Constantan resistance wire, 0.35 mm diameter, 100 m Chrome-nickel resistance wire, 0.25 mm diameter, 100 m Iron resistance wire, 0.2 mm diameter, 100 m Copper (II) sulfate solution 1%, 50 ml	Description         0
Additi Count 1 1 1 1 1 1 1 1 1	CatNo.           520 713           550 42           550 46           550 47           550 51           672 9650           673 5700	uired per class         Name         LIT: LP3 Science Lab Electrics, digital         Constantan resistance wire, 0.35 mm diameter, 100 m         Chrome-nickel resistance wire, 0.25 mm diameter, 100 m         Chrome-nickel resistance wire, 0.35 mm diameter, 100 m         Iron resistance wire, 0.25 mm diameter, 100 m         Copper (II) sulfate solution 1%, 50 ml         Sodium chloride 250 g	Description         Ohmic resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Electrochemistry experiments (LP3.3.7)         Electrochemistry experiments (LP3.3.7)
Additi Count 1 1 1 1 1 1 1 1 1	CatNo. 520 713 550 42 550 46 550 47 550 51 672 9650 673 5700 674 7960	uired per class         Name         LIT: LP3 Science Lab Electrics, digital         Constantan resistance wire, 0.35 mm diameter, 100 m         Chrome-nickel resistance wire, 0.25 mm diameter, 100 m         Chrome-nickel resistance wire, 0.35 mm diameter, 100 m         Iron resistance wire, 0.2 mm diameter, 100 m         Copper (II) sulfate solution 1%, 50 ml         Sodium chloride 250 g         Sulfuric acid, diluted, 0.05 mol/l, 2 l	Description         0         0         0hmic resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Wire's resistance experiments (LP3.3.3)         Electrochemistry experiments (LP3.3.7)         Electrochemistry experiments (LP3.3.7)         Electrochemistry experiments (LP3.3.7)

Detailed information on Mobile-CASSY 2 WiFi, sensors, literature packages and chemical sets are available on the following pages.

### SCIENCE LAB PHYSICS - ELECTRICITY/ELECTRONICS

![](_page_18_Picture_1.jpeg)

With the Mobile-CASSY 2 WiFi, voltage (U), current (I), power (P) and energy (E) can be measured via 4 mm safety sockets.

### LITERATURE PACKAGES

Here you will find an overview of our literature packages. You can find detailed information on our literature on the internet at www.leybold-shop.com.

![](_page_18_Picture_5.jpeg)

![](_page_18_Picture_6.jpeg)

![](_page_18_Picture_7.jpeg)

![](_page_18_Picture_8.jpeg)

# LIT: LP3.3 Electrical basic circuits and electrochemistry

Detailed experiment instructions relating to Science Lab Set EL3 (207 133S). Describes 40 experiments from the field of basic electrical circuits and electrochemistry. <u>Topics:</u> Electrical circuits and switches; Electrical measurement methods; Ohmic resistance; Special resistors; Voltage sources; Electrical application circuits; Electrochemistry

520 7133EN LIT: LP3.3 Electrical basic circuits and electrochemistry

![](_page_18_Picture_12.jpeg)

includes ALL subject areas

### Comprehensive physics experiment instructions in the field of electricity for the Science Lab. Contains 154 experiments on electrostatics, magnetism, basic electrical circuits and electrochemistry, electromagnetism and induction, motors and generators and electronics.

Includes all interactive experiment instructions (Lab Docs) as html file.

![](_page_18_Picture_15.jpeg)

LIT: LP3 Science Lab Electricity, digital

### LIT: LP Science Lab Physics, digital

Comprehensive physics experiment instructions for the Science Lab. Contains 450 experiments in the fields of mechanics, energy, electricity and electronics, optics, atomic and nuclear physics.

Includes all interactive experiment instructions (Lab Docs) as html file.

![](_page_18_Picture_20.jpeg)

### Technical data of the digital version:

- Product key for literature (activation & selection of one literature language in LeyLab)
- Can then be used in LeyLab and Document Center (school/institute licence)
- System requirements:
- Document Center:

- PC with Windows 7 or higher; internet access during installation; local network for distribution to students LeyLab:

- PC, tablet or smartphone with a current browser; internet access

### ADDITIONAL STORAGE ACCESSORIES

![](_page_18_Picture_30.jpeg)

You can find detailed information on additional storage accessories from page 228.

www.ld-didactic.com

### INTRODUCING THE TOPIC

# The possibilities of the plug-in board - changeable and adaptable

![](_page_19_Picture_3.jpeg)

![](_page_19_Picture_4.jpeg)

![](_page_19_Picture_5.jpeg)

FOR SIMPLE ELECTRICAL EXPERIMENTS TO COMPLEX ELECTRONICS CIRCUITS

![](_page_19_Picture_7.jpeg)

### FLEXIBLE AND MODULAR

- Facilitates the use of safety wires in experiments
- Experiment set-up in T and L shapes
- Sturdy
- Minimal space required

![](_page_19_Picture_13.jpeg)

## EASY TO SET UP AND DISASSEMBLE

![](_page_19_Picture_15.jpeg)

The sturdy plug-in boards can be clicked together easily.

## ELECTRICS – EL4

PHYSICS

### OVERVIEW OF OUR CURRICULUM-COMPLIANT EXPERIMENTS

Sensors	LP3.4	ELECTROMAGNETISM AND INDUCTION	
	LP3.4.1	Electromagnetism	
•	LP3.4.1.1 LP3.4.1.2 LP3.4.1.3 LP3.4.1.3C	Magnetic effect of electric current Current-carrying conductor in a magnetic field Magnetic field of a coil Magnetic field of a coil (with Mobile-CASSY 2 WiFi)	DIGITAL
	LP3.4.2	Electromagnetic applications	
	LP3.4.2.1 LP3.4.2.2 LP3.4.2.3	Electromagnet Electromagnetic relays The electric bell	
	LP3.4.3	Induction	
	LP3.4.3.1 LP3.4.3.1C LP3.4.3.2 LP3.4.3.2C	Electromagnetic induction with bar magnet and a coil Electromagnetic induction with bar magnet and a coil (with Mobile-CASSY 2 WiFi) Electromagnetic induction with two coils Electromagnetic induction with two coils (with Mobile-CASSY 2 WiFi)	DIGITAL
	LP3.4.4	Transformers	
•	LP3.4.4.1 LP3.4.4.1C LP3.4.4.2 LP3.4.4.2C	Voltage transformation Voltage transformation (with Mobile-CASSY 2 WiFi) Current transformation Current transformation (with Mobile-CASSY 2 WiFi)	DIGITAL
	LP3.4.5	Applications of induction	
	LP3.4.5.1 LP3.4.5.2 LP3.4.5.2C	Self-induction of a coil (model of an induction coil) Model of an alternating current generator Model of an alternating current generator (with Mobile-CASSY 2 WiFi)	DIGITAL
	LP3.4.6	Coils in direct and alternating current circuits	
	LP3.4.6.1 LP3.4.6.2 LP3.4.6.2C	DC and AC resistance of a coil I (observation experiment) DC and AC resistance of a coil II (measuring experiment) DC and AC resistance of a coil II (measuring experiment) (with Mobile-CASSY 2 WiFi)	DIGITAL
For expe	riments marked wit	th "C", the measurements are carried out digitally with the Mobile-CASSY 2 WiFi.	

🔴 Voltage sensor M, ±30 V 🛛 🔵 Magnetic field sensor M, ±100 mT

![](_page_20_Picture_6.jpeg)

LP3.4.3.2C Electromagnetic induction with two coils

![](_page_21_Picture_1.jpeg)

LP3.4.3.1C Electromagnetic induction with bar magnet and a coil

### OVERVIEW OF EQUIPMENT REQUIRED FOR PERFORMING EXPERIMENTS

![](_page_21_Figure_4.jpeg)

### SCIENCE LAB PHYSICS - ELECTRICITY/ELECTRONICS

![](_page_22_Picture_1.jpeg)

![](_page_22_Picture_2.jpeg)

### Science Lab Electrics EL4 (Set)

Student experiment set of the student experiment system Science Lab in the field of physics. Set-up material for one working group in pre-formed tray. With the equipment set EL4, together with the Science Lab Electrics EL3 (207 133S), 21 experiments at school, college and university level for worldwide curriculums can be performed.

The students deal with the topic electromagnetism. While working out the curriculum required topics, they are also trained in communication and assessment skills. In combination with the Mobile-CASSY 2 WiFi (524 005W), there are additional evaluation options which enable the students digital learning.

### Scope of delivery:

Count	Name	Count	Name
1	Bar magnet	1	Coil 500 turns STE 2/50
1	Plotting compass		Coil 1000 turns STE 2/50
1	Magnetizable rods, set of 4	1	Transformer core, demountable
1	Bell dome	1	Tray, Iow
1	Leaf spring	207 134S	Science Lab Electrics EL4 (Set)

### ADDITIONALLY REQUIRED TO PERFORM ALL EXPERIMENTS

Count	CatNo.	Name	Description	
1	207 133S	Science Lab Electrics EL3 (Set)		
1	521 487	AC/DC Power supply PRO 012 V/3 A		
1	524 005W	Mobile-CASSY 2 WiFi	for digital experiments	
1	524 436	Magnetic field sensor M, ±100 mT		
1	524 438	Voltage sensor M, ±30 V		•
1	500 622	Safety connecting lead 50 cm, blue	Transformation experiment (LP3.4.4)	
2	531 120	Multimeter LDanalog 20	alternative for analog measurements	
Additi	onally requ	uired per class		
Count	CatNo.	Name	Description	
1	520 713	LIT: LP3 Science Lab Electricity, digital		

![](_page_22_Picture_10.jpeg)

# **OVERVIEW OF ADVANTAGES**

- With EL 4, students understand the link between electricity and magnetism for example through induction experiments
- The demountable transformer core is easy to use so the transformation of voltages can be worked on quickly and comprehensibly
- Acquired skills: Understanding the connections between magnetic and electrical phenomena

STUDENT MEASURING	DEVICE	DIGITAL CLASS / EDUCATION
	Mobile-C	ASSY 2 WiFi udent measuring device with WiFi for all measuring tasks in physics, chemistry and biology.
O K	524 005W	Mobile-CASSY 2 WiFi
		You can find detailed information on the Mobile-CASSY 2 WiFi on page 228.
SENSORS		·····
	Magnetic	field sensor M, ±100 mT •

![](_page_23_Picture_7.jpeg)

![](_page_24_Picture_1.jpeg)

### ADDITIONAL STORAGE ACCESSORIES

![](_page_24_Picture_3.jpeg)

You can find detailed information on additional storage accessories from page 228.

### INTRODUCING THE TOPIC

## Generate sounds with electromagnetism

## THE ELECTRIC BELL (LP3.4.2.3)

- Investigating the function of an electric bell
- The materials can be used to set up a bell that is operated with AC voltage

EVERYDAY EXPERIENCE – What happens if I press a bell button? + GAIN KNOWLEDGE ABOUT MAGNETISM + UNDERSTAND ELECTRICAL PHENOMENA

## INTERESTING AND EASY-TO-UNDERSTAND EXPERIMENT INSTRUCTIONS DIGITAL OR IN HARD COPY

Lab Docs are responsive. They adapt the layout to the screen size. From the smallest smartphone tablet or	The electric bell     WORK SHEET       Name:     Date:       Observation     .       0:     Dete:
From the smallest smartphone, tablet or laptop to a projector.	But is Source 3/3 Statistic Source So

# ELECTRICS – EL5

PHYSICS

### OVERVIEW OF OUR CURRICULUM-COMPLIANT EXPERIMENTS

LP3.5	MOTORS AND GENERATORS		
LP3.5.1	Generators		
LP3.5.1.1 LP3.5.1.1C LP3.5.1.2 LP3.5.1.2C LP3.5.1.3 LP3.5.1.3C LP3.5.1.4 LP3.5.1.4C	Dynamo Dynamo (with Mobile-CASSY 2 WiFi) Universal generator - functional principle Universal generator - functional principle (with Mobile-CASSY 2 WiFi) Power plant generator Power plant generator (with Mobile-CASSY 2 WiFi) AC/DC generator with electromagnetic stator AC/DC generator with electromagnetic stator (with Mobile-CASSY 2 WiFi)		DIGITAL DIGITAL DIGITAL
LP3.5.2	Electric motors		
LP3.5.2.1 LP3.5.2.1C LP3.5.2.2	DC motor – functional principle DC motor – functional principle (with Mobile-CASSY 2 WiFi) Universal shunt-wound motor		DIGITAL
LP3.5.2.2C LP3.5.2.3 LP3.5.2.3C	Universal shunt-wound motor (with Mobile-CASSY 2 WiFi) Universal series-wound motor – functional principle Universal series-wound motor – functional principle (with Mobile-CASSY 2 WiFi)	14	DIGITAL
For experiments m	arked with "C", the measurements are carried out digitally with the Mobile-CASSY 2 WiFi.		

![](_page_26_Picture_4.jpeg)

LP3.5.1.1C Dynamo

### SCIENCE LAB PHYSICS - ELECTRICITY/ELECTRONICS

![](_page_27_Picture_1.jpeg)

LP3.5.2.1C DC motor - functional principle

### OVERVIEW OF EQUIPMENT REQUIRED FOR PERFORMING EXPERIMENTS

![](_page_27_Figure_4.jpeg)

### LP3.5 MOTORS AND GENERATORS

### SCIENCE LAB PHYSICS - ELECTRICITY/ELECTRONICS

![](_page_28_Picture_1.jpeg)

![](_page_28_Picture_2.jpeg)

### Science Lab Electrics EL5 (Set)

Student experiment set of the student experiment system Science Lab in the field of physics. Set-up material for one working group. The equipment is stored in Science Lab Electrics EL4 (207 134S). With the supplementary equipment set EL5, together with the Science Lab Electrics EL3 (207 133S) and EL4 (207 134S), 14 experiments at school, college and university level for worldwide curriculums can be performed.

The students deal with the topics motors and generators. While working out the curriculum required topics, they are also trained in communication and assessment skills. In combination with the Mobile-CASSY 2 WiFi (524 005W), there are additional evaluation options which enable the students digital learning.

#### Scope of delivery:

Count	Name	Count	Nan	ne
1	Plug-in axle	1	Coil	rotor STE
1	Rubber rings, set of 8	1 Brush yoke STE		sh yoke STE
1	Pulley Ø 100 mm, plug-in	1	Mag	gneto inductor STE
1	Stator STE 4/50	207 135S		Science Lab Electrics EL5 (Set)

### ADDITIONALLY REQUIRED TO PERFORM ALL EXPERIMENTS

### Additionally required per working group

Count	CatNo.	Name	Description
1	207 133S	Science Lab Electrics EL3 (Set)	
1	207 134S	Science Lab Electrics EL4 (Set)	
1	521 487	AC/DC Power supply PRO 012 V/3 A	
1	524 005W	Mobile-CASSY 2 WiFi	for digital experiments
2	531 120	Multimeter LDanalog 20	alternative for analog measurements
Additio	onally requ	uired per <mark>class</mark>	
Count	CatNo.	Name	Description
	520 713	LIT: LP3 Science Lab Electricity digital	

![](_page_28_Picture_11.jpeg)

# **OVERVIEW OF ADVANTAGES**

- Everyday relevance: Further insights into electromagnetism through experiments with simple motors and generators
- Engine and generator models are quickly assembled and functionally reliable
- Acquired skills: understanding different drive technologies (relevant for the debate on electromobility)

![](_page_29_Figure_5.jpeg)

![](_page_30_Picture_1.jpeg)

![](_page_30_Picture_2.jpeg)

You can find detailed information on additional storage accessories from page 228.

www.ld-didactic.com

### INTRODUCING THE TOPIC

## Electromobility - the e-bike trend

![](_page_31_Picture_3.jpeg)

## CLASSIC: BICYCLE LIGHT WITH DYNAMO

- As with every trend, it's all about gaining basic knowledge
- By spinning the drive roller, the voltage can be measured or visualised using a light bulb
- General understanding of generators
- Transferable to wind energy

![](_page_31_Picture_9.jpeg)

## INNOVATIVE: ELECTRIC DRIVES

- The DC motor is comparable with an e-bike motor
- By assembling and operating different motors, the students gain knowledge of the respective functions
- Through this investigation and further experiments with electric motors, students get their first insight into the topic of electromobility

![](_page_31_Picture_14.jpeg)

# **ELECTRONICS – EL6**

### OVERVIEW OF OUR CURRICULUM-COMPLIANT EXPERIMENTS

LP4.1	BASIC ELECTRONIC CIRCUITS	
LP4.1.1	Capacitors	
LP4.1.1.1 LP4.1.1.1C LP4.1.1.2 LP4.1.1.2C	Capacitors in a DC circuit Capacitors in a DC circuit (with Mobile-CASSY 2 WiFi) Capacitors in an AC circuit Capacitors in an AC circuit (with Mobile-CASSY 2 WiFi)	
LP4.1.2	Relay circuits	
LP4.1.2.1 LP4.1.2.1C	Light-controlled relays Light-controlled relays (with Mobile-CASSY 2 WiFi)	
LP4.1.3	Diodes Characteristic gunua of a dioda	
LP4.1.3.1 LP4.1.3.1C LP4.1.3.2	Characteristic curve of a diode Characteristic curve of a diode (with Mobile-CASSY 2 WiFi) Half-wave rectification	
LP4.1.3.2C	Half-wave rectification (with Mobile-CASSY 2 WiFi)	
LP4.1.3.3C LP4.1.3.4	Full-wave rectification (with Mobile-CASSY 2 WiFi) Light-emitting diodes	
LP4.1.3.4C LP4.1.3.5	Light-emitting diodes (with Mobile-CASSY 2 WiFi) Polarity tester with diodes	
LP4.1.3.5C LP4.1.3.6 LP4.1.3.6	Polarity tester with diodes (with Mobile-CASSY 2 WiFi) Characteristic curve of a Z diode Characteristic curve of a Z diode (with Mobile-CASSY 2 WiFi)	
LP4.1.3.7 LP4.1.3.7C	Overvoltage protection using a Z diode (with Mobile-CASSY 2 Wifi) Overvoltage protection using a Z diode (with Mobile-CASSY 2 WiFi)	
LP4.1.4	Transistors	
LP4.1.4.1 LP4.1.4.2 LP4.1.4.2C	Diode junctions on transistors, test circuit with light-emitting diodes Transfer characteristic of a transistor Transfer characteristic of a transistor (with Mobile-CASSY 2 WiFi)	
LP4.1.4.3 LP4.1.4.3C LP4.1.4.4	Transistor circuit I, voltage control Transistor circuit I, voltage control (with Mobile-CASSY 2 WiFi) Light-controlled transistor I, light barrier	
LP4.1.4.5 LP4.1.4.5C	Delay switch Delay switch (with Mobile-CASSY 2 WiFi)	
LP4.1.5	Diode circuits	
LP4.1.5.1 LP4.1.5.1C	Overvoltage and reverse polarity protection using diodes Overvoltage and reverse polarity protection using diodes (with Mobile-CASSY 2 WiFi)	
LP4.1.5.2	Smoothing pulsating DC voltages with capacitors	
LP4.1.6	Flip-flops	
LP4.1.6.1 LP4.1.6.1C	Basic experiments with flip-flops Basic experiments with flip-flops (with Mobile-CASSY 2 WiFi)	
LP4.1.7	Amplifier circuits	
LP4.1.7.2 LP4.1.7.2C	Touch-sensitive switches, humidity and fill level indicators Touch-sensitive switches, humidity and fill level indicators (with Mobile-CASSY 2 WiFi)	
LP4.1.8	Solar cells	
LP4.1.8.1 LP4.1.8.1C	Forward and reverse direction of a solar cell Forward and reverse direction of a solar cell (with Mobile-CASSY 2 WiFi)	
LP4.1.8.2	Output and power characteristics of a solar cell	
LP4.1.8.2C	No-load voltage of a solar cell	
LI T. I.O.J		

www.ld-didactic.com

![](_page_33_Picture_1.jpeg)

```
LP4.1.3.3C Full-wave rectification
```

### OVERVIEW OF EQUIPMENT REQUIRED FOR PERFORMING EXPERIMENTS

![](_page_33_Figure_4.jpeg)

### SCIENCE LAB PHYSICS - ELECTRICITY/ELECTRONICS

![](_page_34_Picture_1.jpeg)

### Science Lab Electrics EL6 (Set)

Student experiment set of the student experiment system Science Lab in the field of physics. Set-up material for one working group in pre-formed tray. With the equipment set EL6, together with the Science Lab Electrics EL3 (207 133S), 42 experiments at school, college and university level for worldwide curriculums can be performed.

The students deal with the topics basic electronic circuits and transistor electronics. While working out the curriculum required topics, they are also trained in communication and assessment skills. In combination with the Mobile-CASSY 2 WiFi (524 005W), there are additional evaluation options which enable the students digital learning.

#### Scope of delivery:

Count	Name
1	Plug-in board safety socket, 20/10
1	Resistor 470 Ohm, STE 2/19
1	Resistor 4.7 kOhm, STE 2/19
1	Potentiometer 220 Ohm, STE 4/50
1	Capacitor, 1 µF, STE 2/19
1	Capacitor (electrolytic) 100 $\mu$ F, STE 2/19
1	Capacitor (electrolytic) 470 µF, STE 2/19
1	Light emitting diode red, STE 2/19
4	Diode 1N 4007, STE 2/19
1	Zener diode 6.2, STE 2/19

Count	Name		
1	Light emitting diode green, STE 2/19		
1	Photodiode, lateral		
1	Solar cell, STE 2/19		
1	Transistor BD 137, NPN, e.b., STE 4/50		
1	Transistor BD 138, PNP, e.b., STE 4/50		
1	Relay with change-over switch STE 4/50		
1	Earphone		
1	Tray, low		
207 1365	Science Lab Electrics EL6 (Set)		

### ADDITIONALLY REQUIRED TO PERFORM ALL EXPERIMENTS

Count	CatNo.	Name	Description		
1	207 133S	Science Lab Electrics EL3 (Set)			
1	521 487	AC/DC Power supply PRO 012 V/3 A			
1	524 005W	Mobile-CASSY 2 WiFi	for digital experiments		
1	524 438	Voltage sensor M, ±30 V	•		
1	500 622	Safety connecting lead 50 cm, blue	Overvoltage protection experiment (LP4.1.3)		
2	531 120	Multimeter LDanalog 20	alternative for analog measurements		
Additionally required per class					
Count	CatNo.	Name	Description		
1	520 713	LIT: LP3 Science Lab Electricity, digital			

![](_page_34_Picture_10.jpeg)

PHYSICS

# **OVERVIEW OF ADVANTAGES**

- Students learn to understand more complex structures, such as diode and transistor circuits, through the structured experiment instructions
- Similarities to the use of printed circuit boards in electronic components become visible
- Acquired skills: understanding the functions of electronic components in modern technical devices

STUDENT MEASURING	DEVICE	DIGITAL CLASS / EDUCATION				
	Mobile-(	CASSY 2 WiFi				
A THE A	The universal s	tudent measuring device with WiFi for all measuring tasks in physics, chemistry and biology.				
O X	524 005W	Mobile-CASSY 2 WiFi				
		You can find detailed information on the Mobile-CASSY 2 WiFi on page 228.				
SENSORS						
	For measuring integrated vol	the electrical voltage up to $\pm 30$ V with Mobile-CASSY 2 WiFi (524 005W). In connection with the tage input, Mobile-CASSY 2 WiFi (524 005W) can become a 2 channel storage oscilloscope.				
	524 438	Voltage sensor M, ±30 V				
		You can find detailed information on this and other sensors from page 229.				

![](_page_36_Picture_1.jpeg)

### ADDITIONAL STORAGE ACCESSORIES

![](_page_36_Picture_3.jpeg)

You can find detailed information on additional storage accessories from page 228.

### INTRODUCING THE TOPIC

# Ohm's Law as a diagram with the Mobile-CASSY 2 WiFi

![](_page_37_Figure_3.jpeg)

#### Current and voltage measurement on a resistor

## MEASURE SEVERAL PHYSICAL VARIABLES AT THE SAME TIME

- With the Mobile-CASSY 2 WiFi, two or more measuring quantities can be measured against each other and simultaneously recorded, such as:
  - Current and voltage
  - Temperature and voltage
- Dependencies between measuring quantities can be illustrated exceptionally well in diagrams
- This allows students to gain a deeper understanding of the basic principles of electronics
- Simple measuring of characteristic curves is only possible with digital measuring technology

## DIRECT DISPLAY OF CHARACTERISTIC CURVES IN THE DIAGRAMS IN THE DIGITAL EXPERIMENT LITERATURE

With interactivity between the Mobile-CASSY 2 WiFi and Lab Docs, the measured values are transferred in real time and the characteristic curves are displayed directly in the diagrams.

The Lab Doc with the recorded characteristic curve can also be saved as a digital protocol and then shared with the teacher.

With the Lab Docs Editor, the diagrams (among other features) can be adjusted.

![](_page_37_Picture_16.jpeg)